



You can find out within minutes whether or not you are expecting a baby, but as you await this life-changing result, it can seem like a lifetime before the symbol appears.

HOW TO USE A HOME PREGNANCY TEST

Always read the instructions, but most tests work as follows:

- **You urinate on the stick** and leave it for a specified number of minutes.
- **A symbol will appear** in the control window to indicate the test is working (if this does not appear, the test is faulty). If a symbol then appears in the results window, you are pregnant.
- **It is advisable to do the test** first thing in the morning when your urine is most concentrated. By doing this, the hCG levels are likely to be detected more easily.
- **The results symbol** gradually fades, so read it after the specified time, not later. Do another test the next morning if you're unsure of the result.